

# Varnal<sup>®</sup>

## CAVITY VARNISH

Varnal is a superior cavity varnish and tubuli seal. Varnal dries fast, hard, and is visible as a dull frosting on tooth surfaces when application is completed. Easily applied with a cotton pledget or brush, Varnal requires only one application. Varnal's special solvent system retains solids in solution; therefore no separate thinner is ever required. Varnal ensures restoration success and tooth protection from amalgams and other restoration materials.



Description	Size	Item No.
Varnal <sup>®</sup> Cavity Varnish	14 g	0801

## Key Features

- One-step fast-drying cavity varnish
- Protects dentin, penetrates, fills, and seals tubuli
- Only one application required per tooth preparation
- Tooth need not be dried prior to application
- No burnishing required
- Compatible with amalgams and other restoration materials

## Benefits

- Won't interfere with patient eating or drinking after procedure
- No thinner required to maintain original viscosity

Cetylite Inc.  
9051 River Road  
Pennsauken, NJ 08110  
856-665-6111  
marketing@cetylite.com

©2017 Cetylite Industries, Inc. All rights reserved. Information is summary in nature and subject to change. Varnal and Cetylite are registered trademarks of Cetylite Industries, Inc.

01CETY2017